

12-14-04 3P5H10 IFL on Human Tissues for Montse/Steve Finebeiner
(Paraffin Embedded)

1:3000 }
1:1000 } 3 Control Cerebellums + 3 affected (blinded)
1:300 } = six slides

also 1:5000 CB Rb as a control (A) SE6 SE18
0.1M glycine - 0.3M

1. 60°C to melt paraffin - 20 mins
2. deparaffinize 3x xylenes 5 min each in 1M Tris base
100% EtOH " 15-30 mins, wash
95% EtOH " w/PBS
70% EtOH "

3. Tritox x-100 0.5% in PBS 2 x 15-20 min ; Tris wash - 5 min

4. Glycine Block - 0.3M in Tris 0.1M - 45-60 min
75.07g x 0.3 mol x 1L x 30 ml = 0.67563 g in 30 ml Tris 0.1M
mol 1000ml

5. rinse w/ PBS - 10 min, PBS-Tx - 10 min

6. block - NDS - 10% - 2.5 ml .. 60 min - 120 min
25mls gel - 0.2% - 500 ul

milk - 2% - 500 mg
PBS-Tx - 22 ml 0.5% Triton in PBS-Tx 5ml

7. rinse PBS-Tx 10 min

8. 1° incubation O/N RT 1:300 - 3ul in 300ul 150ul - 3% NDS
150ul 1:1000 - 30ul (1:100) in 270ul 0.2% of lectin
CB - 1:5000 1:3000 - 10ul (1:100) in 290ul PBS-Tx 4.75ml
6 ul (1:100) for 300 ul

9. rinse - PBSTx 15 min, PBS 15 min, PBS-Tx 15 min, PBS 15 min, PBSTx 15 min

10. 2° incubation 1:300 cy2 x Ms (DK) - 3P5H10 - 10ul

diluent: 1:300 cy3 x Rb (DK) - Calbindin - 10ul

3% NDS - 90ul

gelatin - 60ul

PBS-Tx - 2.83 ml

11. rinses

(A) if not fragile tissues

1% Na Borohyd. in PBS for 30 min

wash w/ PBS 15 min
5 min

RESULTS. possibly too much auto fluorescence

glycine block not rigorous enough

Next: try 1% sodium borohydride / PBS

proper controls.

No 1° or 2°

No 1°

Both 1° + 2°

10 mM citrate buffer pH 6.0 15 min

90% formic acid 7 min

SIGNATURE

READ AND UNDERSTOOD

DATE 20

DATE 20

Serum

dk x Ms
goat x Rb

wash 1x PBS 10 min each

3x 10 min